

AUT Tables

One table, multiple solutions



Fast & Professional

About our AUT Tables

Lab Design by Fenco Solution's Adjustable Utility Tables allow you the flexibility to design your workspace exactly how you want it. From electrical and plumbing fixtures to height adjustability and storage, look no further than our AUTs. Our heavy-duty table systems let you lay out your laboratory incorporating versatile arrangements that can quickly adjust to your business' evolving needs.









www.FENCOLABS.com

Find your solution today!



Toll-Free: 800.486.8484

Features & Options

-  **Adjustable height**
(30"-36" in 1" increments)
-  **Stationary or mobile options**
-  **Up to 3 plumbing fixtures**
pre-plumbed on each upright
-  **3 pre-wired duplex power strip**
Additional duplex can be added upon request
-  **14 GA. painted metal shelves**
Additional shelving material options are available
-  **Available in all 20+ Standard Fenco Colors**



Sizes

- Widths: 36", 48", 60", 72"
- Heights: 78"
- Depth: 30"

Lead Time

- 8-10 Weeks (Fully Customized)

Storage

Mobile or suspended casework brings additional storage to your AUT. Adjustable shelving with seismic lips ensures your material stays safely stored above your counter.

Design

AUTs were engineered with the user in mind. Designed and fabricated in the US, this workstation offers a range of ergonomic and flexible elements to ensure your workstation suits your requirements.

Additional Technical Writeup

AUTs are available in four common widths - 36", 48", 60" and 72" - with a standard height of 78"H. Enhance the mobility of your table by incorporating casters, allowing easy movement within your lab. Each table can accommodate a maximum of six pre-plumbed service fixtures. The UL listed electrical housings come standard, with 6 outlets and 10' long power cord. You have the flexibility to choose a NEMA plug and outlet type according to your specific needs.

Ordering made easy

- 1 Determine cabinet material
- 2 Select measurements of cabinet
- 3 Specify drawer/door configurations
- 4 Select Color



www.FENCOLABS.com

Find your solution today!



Toll-Free: 800.486.8484


ADJUSTABLE UTILITY TABLES

AUT-4830S, 6030S, 7230S


Ordering Information

- For Casters, add suffix 'C' to end of model number (Ex. 4830C)
- For Electrical Power Strip with (3) duplexes, add suffix '-E' to end of model number (Ex. 4830C-E)
- Add suspended or mobile casework to meet your storage needs.
- Electrical power strip comes standard with 20A twist lock plug.

AUT-6030 with Casters



AUT-6030



AUT Tables w/ no Gases	AUT Tables w/ 1 Gas (G1)	AUT Tables w/ 2 Gases (G2)	AUT Tables w/ 3 Gases (G3)
	Duplex (3)	Duplex (3)	Duplex (3)
	L4200-651 Ball Valve and Hose with Quick Connect (1)	L4200-651 Ball Valve and Hose with Quick Connect (2)	L4200-651 Ball Valve and Hose with Quick Connect (2)
	Pre-plumbed (specify gas)	Pre-plumbed (specify gas)	Pre-plumbed (specify gas)
48"W 30"D 12"D Shelves (2)			
AUT-4830.....	AUT-4830-G1.....	AUT-4830-G2.....	AUT-4830-G3.....
AUT-4830C.....	AUT-4830C-G1.....	AUT-4830C-G2.....	AUT-4830C-G3.....
AUT-4830C-E.....	AUT-4830C-E-G1.....	AUT-4830C-E-G2.....	AUT-4830C-E-G3.....
AUT-4830-E.....	AUT-4830-E-G1.....	AUT-4830-E-G2.....	AUT-4830-E-G3.....
60"W 30"D 12"D Shelves (4)			
AUT-6030.....	AUT-6030-G1.....	AUT-6030-G2.....	AUT-6030-G3.....
AUT-6030C.....	AUT-6030C-G1.....	AUT-6030C-G2.....	AUT-6030C-G3.....
AUT-6030C-E.....	AUT-6030C-E-G1.....	AUT-6030C-E-G2.....	AUT-6030C-E-G3.....
AUT-6030-E.....	AUT-6030-E-G1.....	AUT-6030-E-G2.....	AUT-6030-E-G3.....
72"W 30"D 12"D Shelves (4)			
AUT-7230.....	AUT-7230-G1.....	AUT-7230-G2.....	AUT-7230-G3.....
AUT-7230C.....	AUT-7230C-G1.....	AUT-7230C-G2.....	AUT-7230C-G3.....
AUT-7230C-E.....	AUT-7230C-E-G1.....	AUT-7230C-E-G2.....	AUT-7230C-E-G3.....
AUT-7230-E.....	AUT-7230-E-G1.....	AUT-7230-E-G2.....	AUT-7230-E-G3.....